

Automated Detection and Localization of Abnormalities in Chest X-rays using Weakly-Supervised Learning

Presentation order: 15

Members:

Lu Jinzhou A0329981M

Song Yiran A0294421H

Wu Yuhang A0294425Y

Introduction

- **Chest X-rays** are essential for diagnosing conditions like pneumonia, pleural effusion, and cardiomegaly.
- Traditional interpretation is **time-consuming** and relies heavily on **physician experience**.



Cord Challenge:

- **High Cost of Data Annotation:** Pixel-level annotation for medical images is expensive and time-consuming, a major bottleneck for AI applications.
- "Black Box" Model Problem: The opaque decision-making process of deep learning models leads to a lack of trust from clinicians.

Goal:

 To develop an AI-assisted diagnostic framework with both high accuracy and high interpretability to address the core challenges of data cost and model trust.

Problem Statement

Task

Multi-Label Classification & Explainable AI (XAI)
 Develop a model that not only predicts multiple chest abnormalities but also explains its decisions.

Input & Process

- Input: A single frontal chest X-ray image from the MIMIC-CXR dataset.
- **Process:** Train a deep convolutional neural network using only image-level "weak labels" (e.g., "Pneumonia").

Output

- Classification: Predict the presence of multiple common thoracic abnormalities.
- Localization & Explanation: Generate a heatmap to visually highlight regions the model focuses on for diagnosis.

Methodology

Data Processing

- **Dataset:** Subset of the public MIMIC-CXR dataset (1000 samples), which contains chest X-rays and their corresponding radiology reports. Split into training, validation and test sets with a 8:1:1 ratio.
- Automated Weak Labeling: Use the CheXpert NLP tool to automatically extract image-level labels from free-text radiology reports.
- **Advantage:** Bypasses the need for manual pixel-level annotation, significantly improving data preparation efficiency.
- **Preprocessing:** Normalize all images to a fixed size (e.g., 224x224 pixels), transformed into tensors and normalized using ImageNet's mean and standard deviation values

Model Selection

- **Backbone Model:** primarily use DenseNet-121 with ResNet-50 as comparison, pre-trained on ImageNet.
- **Reasoning:** Achieves high performance with greater parameter efficiency, which helps reduce the risk of overfitting.


Methodology-Comparison


DenseNet-121 Architecture Features

- Dense Connectivity: Each layer is connected to every other layer in a feed-forward fashion, which strengthens feature propagation and encourages feature reuse.
- Parameter Efficiency: Achieves high performance with fewer parameters than traditional CNNs, reducing the risk of overfitting.
- Mitigates Vanishing Gradients: The dense connection path improves the flow of information and gradients throughout the network.

ResNet-50 Architecture Features

- Residual Learning: Introduces "shortcut connections" that allow the network to learn residual functions, making it easier to train much deeper networks.
- Solves Degradation Problem: Effectively addresses the performance degradation issue in deep networks, allowing for a significant increase in the number of layers.
- Modular Structure: Composed of stacked residual blocks, making the architecture clean and easy to extend.






Methodology

Training Methodology

- Loss function: BCEWithLogitsLoss (Suitable for multi-label classification).
- Optimizer: Adam Optimizer.
- Output Layer: A Sigmoid activation function is used to support multi-label classification (i.e., one image can have multiple findings).
- XAI(Explainable AI): Apply Grad-CAM (Gradient-weighted Class Activation Mapping)
 post-training to generate explanation heatmaps.
- Auxiliary Methods: Explore methods like LIME, DeepLift, and Integrated Gradients to understand the model's behavior from different perspectives and enhance trust in the results.

Evaluation Metrics


• Quantitative Evaluation: We use standard metrics such as AUROC, Accuracy, F1-Score, and confusion matrices on an independent test set to evaluate classification performance.



Initial Progress

Preliminary Data Analysis

- Data pre-processing: extract DICOM metadata, organize anonymized reports, and export CSV files;
- **NLP report labeling:** applied the CheXpert pipeline to classify 14 chest findings;
- **Statistical analysis:** generate descriptive statistics and label distributions;
- Visualization & reporting



Sop 10 - PeoPrive Calest Court

934
Total Reports Processed

242Total Unique Patients

2732Total Positive Findings

Summary & Novelty

Summary:

Our project provides a practical solution for developing efficient, trustworthy,
 and implementable clinical AI tools.

Novelty:

- Automated Weak Labeling: Uses CheXpert NLP to generate multi-label image annotations directly from radiology reports, bypassing manual labeling.
- **Efficient Multi-label Classification**: DenseNet-121 backbone with Sigmoid outputs enables accurate prediction of multiple findings per X-ray.
- **Explainable AI Integration**: Grad-CAM and auxiliary XAI methods provide visual explanations, improving model interpretability and trust.

Reference:

[1]Tanno, R., Barrett, D. G. T., 2025. Collaboration between clinicians and vision-language models in radiology report generation. Nature medicine, 31(2), 599–608. https://doi.org/10.1038/s41591-024-03302-1
[2]L. N. Rohmah and A. Bustamam, "Improved Classification of Coronavirus Disease (COVID-19) based on Combination of Texture Features using CT Scan and X-ray Images," 2020 3rd International Conference on Information and Communications Technology (ICOIACT), Yogyakarta, Indonesia, 2020, pp. 105-109, doi: 10.1109/ICOIACT50329.2020.9332123.
[3]Thapa, M., & Kaur, R. (2025). An explainable deep-learning based multi-label image classification for chest X-rays. Procedia Computer Science, 239, 281–288.https://doi.org/10.1016/j.procs.2025.04.505

Code Lab(Up to date):

https://github.com/Yuhang20/BN5212.git https://github.com/stanfordmlgroup/chexpert-labeler.git

THANK YOU